



VAH listed

Surface Des

Concentrated disinfectant and cleaner for surfaces

- user-friendly active ingredients
- strong cleaning properties
- suitable for acrylic glass
- perfume-free



Product properties

Adelhyde-free concentrate for disinfecting and cleaning of medical devices, medical inventory and surfaces of all sorts. The product is perfume-free and consequently appropriate for use in food product industries.

Composition

100 g solution contain:
 6,66 g alkyl(C12-16)dimethylbenzyl ammonium chloride
 3,33 g Dimethyldidecyl ammonium chloride

Chemical physical data

Appearance of concentrate: clear, colourless liquid
 PH-value: 7,25 +/- 1,25

Range of Application

Acc. to EU Medical Devices Regulation:
 Rapid disinfection and cleaning of non-invasive alcohol-sensitive medical devices.
Acc. to Biocidal Products Regulation:
 Disinfection and cleaning of all kinds of surfaces and medical inventory. Also for all areas in the food industry and canteen kitchens.

Application

Surface wipe disinfection: Fully wet the surfaces to be disinfected with a sufficient quantity of the prepared working solution and spread mechanically (scrub-wipe method).

Concentration and exposure times

Application recommendation for surface disinfection:	
bactericidal, yeasticidal:	7.5 % (75 g/L) – 5 min 2.5 % (25 g/L) – 15 min 1 % (10 g/L) – 60 min
active against all enveloped viruses:	2.5 % (25 g/L) – 5 min 1 % (10 g/L) – 60 min
Papova-/Polyomavirus:	7.5 % (75 g/L) – 2 min 2 % (20 g/L) – 30 min
rotavirus:	0.5 % (5 g/L) – 30 s

Delivery Forms

Single package 5 Litre canister	Sales unit 1 x 5 L	Ref.-No. 09941
------------------------------------	-----------------------	-------------------

Pharmaceutical expert report

Prof. Dr. Schubert, Frankfurt/Main:
 Certificate for Surface Disinfection according to VAH Guidelines from 15.04.2010

Prof. Dr. Werner, Bischofshofen:
 Certificate for Surface Disinfection according to VAH Guidelines from 13.04.2010

Dr. Steinmann, Virologe, Bremen:
 Certificate for Vaccinia-virus effectiveness from 13.08.2010,
 Certificate for BVDV effectiveness from 24.08.2010,
 Certificate for Rota-virus effectiveness from 12.12.2010

Dr. Suchomel, Wien:
 Certificate for Surface Disinfection according to VAH Guidelines from 28.02.2017

Classification according regulation (EC) Nr. 1272/2008:



Danger: Quaternary ammonium compounds, benzyl-C12-16-alkyldimethyl, chlorides; didecyldimethylammonium chloride; alcohols, C9-11-branched, ethoxylated.
 Causes skin irritation. Causes serious eye damage. Very toxic to aquatic life with long lasting effects. Wear protective gloves/ protective clothing/eye protection/face protection. IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Avoid release to the environment. Immediately call a POISON CENTER/ doctor. Dispose of contents/container to approved disposal company or local collection. For professional use only.

Use disinfectants in a safe way. Always read label and product information before use.

CE 0482 UN 1903
 baua: Reg.-Nr. N-50653, N-57130



AMPri Handelsgesellschaft mbH
 Benzstr. 16 / 21423 Winsen (Luhe), Germany
 Tel. +49 (0)4171 - 8480-0 / www.ampri.de / info@ampri.de

bioSight GmbH
 Brasselsbergstr. 3 / 34132 Kassel, Germany
 Tel.: +49 (0)561 - 81678111