

## Technical Data Sheet

Article-No.: **05563-Y**  
Description: **MED-COMFORT**  
non woven protective coat  
with PE-coating



SIZES	
	L/XL
width	138 cm
length	138 cm

PRODUCT DESCRIPTION							
material	<input checked="" type="checkbox"/> Polypropylene (PP)	<input checked="" type="checkbox"/> Polyethylene (PE)	<input type="checkbox"/> SMS				
colour	<input type="checkbox"/> blue [B]	<input type="checkbox"/> green [G]	<input type="checkbox"/> orange [O]	<input type="checkbox"/> red [R]	<input type="checkbox"/> white [W]	<input checked="" type="checkbox"/> yellow [Y]	<input type="checkbox"/> light-blue [LB]
	<input type="checkbox"/> dark-blue [DB]	<input type="checkbox"/> dark-green [DG]					
closure	<input checked="" type="checkbox"/> to tie	<input type="checkbox"/> zipper	<input type="checkbox"/> push buttons	<input type="checkbox"/> velcro			
characteristics	<input checked="" type="checkbox"/> tricot-cuff	<input type="checkbox"/> elastic wrist-cuff	<input checked="" type="checkbox"/> latex-free	<input type="checkbox"/> hood	<input type="checkbox"/> stand-up collar	<input checked="" type="checkbox"/> partially laminated	<input type="checkbox"/> sewed PE-oversleeves
	<input type="checkbox"/> shortenable by perforation	<input type="checkbox"/> blocked	<input type="checkbox"/> dust-cover	<input type="checkbox"/> embossed	<input type="checkbox"/> short sleeves	<input type="checkbox"/> with cardboard	<input checked="" type="checkbox"/> single folded

REGULATORY AFFAIRS							
PPE-Regulation (EU) 2016/425	<input type="checkbox"/> category I	<input type="checkbox"/> category II	<input checked="" type="checkbox"/> category III	<input type="checkbox"/> no PPE-article			CE
MD-Regulation (EU) 2017/745	<input type="checkbox"/> class I	<input type="checkbox"/> class II	<input type="checkbox"/> class III	<input type="checkbox"/> sterile	<input type="checkbox"/> measuring function	<input checked="" type="checkbox"/> no medical device	
Food Contact (EG) 1935/2004	<input checked="" type="checkbox"/> acidic foods	<input checked="" type="checkbox"/> aqueous foods	<input checked="" type="checkbox"/> fatty foods	<input checked="" type="checkbox"/> alcoholic foods	<input checked="" type="checkbox"/> dry foods	<input type="checkbox"/> not approved for food-contact	

STANDARDS							
tested in accordance to	<input checked="" type="checkbox"/> DIN EN 14126	<input checked="" type="checkbox"/> ISO 16603	<input checked="" type="checkbox"/> ISO 16604	<input checked="" type="checkbox"/> ISO 22610	<input checked="" type="checkbox"/> ISO 22611	<input checked="" type="checkbox"/> ISO 22612	
	<input checked="" type="checkbox"/> EN ISO 13688	<input checked="" type="checkbox"/> EN 13034					

## Technical Data Sheet

Article-No.: **05563-Y**  
Description: **MED-COMFORT**  
non woven protective coat  
with PE-coating



LOGISTIC DATA SUBPACKING	
general information	
material	bag
pieces per subpacking	10
GTIN subpacking blue [B]	-
GTIN subpacking green [G]	-
GTIN subpacking orange [O]	-
GTIN subpacking red [R]	-
GTIN subpacking white [W]	-
GTIN subpacking yellow [Y]	4044941753560
GTIN subpacking light-blue [LB]	-
GTIN subpacking dark-blue [DB]	-
GTIN subpacking dark-green [DG]	-
measures & sizes	
length	430 mm
width	330 mm
height	-
weights	
size	gross weight*
L/XL	985 g

LOGISTIC DATA PALETTE	
general information	
kind of palette	euro-palette
measures & sizes	
cartons per layer	18
layer per palette	#BEZUG!
palette height	174
weights	
size	gross weight*
L/XL	211 kg

LOGISTIC DATA OUTER PACKING	
general information	
material	carton
subpackings per outer packing	10
GTIN outer packing blue [B]	-
GTIN outer packing green [G]	-
GTIN outer packing orange [O]	-
GTIN outer packing red [R]	-
GTIN outer packing white [W]	-
GTIN outer packing yellow [Y]	4044941753577
GTIN outer packing light-blue [LB]	-
GTIN outer packing dark-blue [DB]	-
GTIN outer packing dark-green [DG]	-
measures & sizes	
length	600 mm
width	470 mm
height	265 mm
weights	
size	gross weight*
L/XL	10.350 g

## Technical Data Sheet

Article-No.: **05563-Y**  
Description: **MED-COMFORT**  
**non woven protective coat**  
**with PE-coating**



### WARNINGS AND SAFETY INFORMATION

<b>storage / expiry date</b>	Store product in original packaging in a cool and dry place without additional weight, protect from direct sunlight. Do not store near ozone sources (laser printers, copiers). The actual expiry time in use cannot be specified in general terms, as it depends on the general conditions of use. An individual risk assessment must be carried out in each case. The expiry date - valid for proper storage - is stated on the packaging.
<b>use and control</b>	Only ever use the product for its intended use. A check/risk assessment must be carried out to ensure that the product is suitable for the intended use. Ensure that the size matches that of the user. Do not make any changes to the product. Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.). This product is intended for single use only
<b>disposal</b>	The used product must be disposed of in accordance with the disposal regulations of the local waste disposal company. Unused products can be disposed of with household waste.
<b>warnings</b>	Choose products that are compatible with the work area. - The disposable item must be replaced after each use. - In the event of tears, punctures, etc., leave the work area and wear a new protective gown. - Prolonged wearing of protective gowns can lead to heat stress. Heat stress and discomfort can be reduced or eliminated by using suitable underwear or suitable ventilation equipment. - These garments are flammable - keep away from fire. - In the event of damage to the product, leave the workplace immediately. - The user must not remove the protective gown while still in the danger zone.
<b>donning and doffing instructions</b>	<p><b>WAY OF PUTTING ON:</b></p> <ul style="list-style-type: none"> <li>- Make sure that the size matches that of the user. Do not make any changes to the product.</li> <li>- Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.).</li> <li>- Open the gown and put it on. Make sure that the garment does not tear. Close the gown by wrapping the waistband around your body and fastening it with a bow.</li> </ul> <p><b>REMOVAL METHOD:</b></p> <ul style="list-style-type: none"> <li>- Undo the waist and collar ties and remove the gown.</li> <li>- The disposal of contaminated clothing is regulated by national or local laws.</li> </ul> <p>Do not dispose of contaminated protective clothing in household waste.</p> <ul style="list-style-type: none"> <li>- The protective properties are only effective if the article is correctly put on and closed.</li> <li>- Protect uncovered parts of the body (hands, respiratory tract, feet) with protective gloves, boots, mask, etc., attach them to the gown (with adhesive strips if necessary) and ensure that they offer the same level of protection to ensure full-body protection.</li> </ul>

\*slight deviations possible due to standard tolerances

rev-no.: 2025-01  
date: 26.06.2025

changes and errors excepted