

### AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

# **Technical Data Sheet**

Article-No.: 02420

Description: STYLE by Med-Comfort

facemask

14683:2019

accordance to


EN

3-ply, with elastic-band



length	width						
175 mm	95 mm						
PRODUCT D	DESCRIPTION	_					
material	✓ non-woven PP	□ paper					
colour	☐ blue [B]	☐ berry [Berry]	☐ bordeaux [BR]	cedro [C]	green [G]	☐ lagoon [Lagoon]	purple [L]
	☐ lime [Lime]	☐ Mix (blue,yellow,	☐ orange [O]	☐ pink [P]	☐ white [W]	yellow [Y]	☐ platinum / grey [PT]
	<ul><li>✓ blue patterned</li><li>[DB]</li></ul>	wine-red patterned [WR]	☑ black [SCH]	☐ black-red- gold	clean ocean / ocean green [CO]		
characteristic	type II	☑ 3 layers	☑ latex-free	☐ with valve	☐ anti-fog	with ear loops	☐ with activated carbon filter
characteristic	s bacterial filtration	on capacity ≥ 98%					
REGULATO	RY AFFAIRS						
PPE-Regulation (EU) 2016/42	on	☐ category II	☐ category III	ono PPE- article			CE
MD-Regulatio (EU) 2017/74		☐ class II	☐ class III	☐ sterile	<ul><li>measuring fundtion</li></ul>	no medical device	00
STANDARD	<u></u>						
tested in	□ EN	EN 149:2001					

QMFORM\_67.001 1/3 Herausgabe: 15.09.2021





Benzstraße 16 21423 Winsen (Luhe) - Germany

> Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

# **Technical Data Sheet**

Article-No.: **02420** 

Description: STYLE by Med-Comfort

facemask

3-ply, with elastic-band



LOGISTIC DATA SUBPACKING	
general information	
material	carton
pieces per subpacking	50
GTIN subpacking blue [B]	-
GTIN subpacking berry [B]	-
GTIN subpacking bordeaux [BR]	-
GTIN subpacking cedro [C]	-
GTIN subpacking green [G]	-
GTIN subpacking lagoon [Lagoon]	-
GTIN subpacking purple [L]	-
GTIN subpacking lime [L]	-
GTIN subpacking mix [Mix]	-
GTIN subpacking orange [O]	-
GTIN subpacking pink [P]	-
GTIN subpacking white [W]	-
GTIN subpacking yellow [Y]	-
GTIN subpacking platinum [PT]	-
GTIN subpacking blue patterned [DB]	4044941031941
GTIN subpacking wine-red patterned [WR]	4044941031385
GTIN subpacking black [SCH]	-
GTIN subpacking black-red-gold	-
GTIN subpacking clean ocean [CO]	-
measures & sizes	
length	180 mm
width	100 mm
height	90 mm
weights	
size	brut weight*
uni	165 g

LOGISTIC DATA PALETTE	
general information	
kind of palette	euro-palette
measures & sizes	
cartons per layer	8
layer per palette	8
palette height	175
weights	·
size	brut weight*
uni	163 kg

LOGISTIC DATA OUTER PACKING	
general information	
material	carton
subpackings per outer packing	10
GTIN outer packing blue [B]	-
GTIN outer packing berry [B]	-
GTIN outer packing bordeaux [BR]	-
GTIN outer packing cedro [C]	-
GTIN outer packing green [G]	-
GTIN outer packing lagoon [Lagoon]	-
GTIN outer packing purple [L]	-
GTIN outer packing lime [L]	-
GTIN outer packing mix [Mix]	-
GTIN outer packing orange [O]	-
GTIN outer packing pink [P]	-
GTIN outer packing white [W]	-
GTIN outer packing yellow [Y]	-
GTIN outer packing platinum [PT]	-
GTIN outer packing blue patterned [DB]	4044941031958
GTIN outer packing wine-red patterned [WF	4044941031392
GTIN outer packing black [SCH]	-
GTIN outer packing black-red-gold	-
GTIN outer packing clean ocean [CO]	-
measures & sizes	
length	520 mm
width	190 mm
height	200 mm
weights	
size	brut weight*
uni	2.150 g



### AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

## **Technical Data Sheet**

Article-No.: 02420

Description: STYLE by Med-Comfort

facemask

3-ply, with elastic-band





### WARNINGS AND SAFETY INFORMATION

storage / expir	Store face mask in original packaging in a cool and dry place without additional weight, protect from direct sunlight. Do
date	not store near ozone sources (laser printers, copiers). The actual expiry time in use cannot be specified in general terms,
	as it depends on the general conditions of use. An individual risk assessment must be carried out in each case. The
	expiry date - valid for proper storage - is stated on the packaging.

use and control	Always use face masks only for their intended use. A check/risk assessment must be carried out to ensure that the face
	mask is suitable for the intended use.
	Change the mask when it is soaked through.

disposal	Used face masks must be disposed of in accordance with the disposal regulations of the local waste disposal company.
	Unused face masks can be disposed of with household waste.

warnings	Change the mask when it is soaked through.  Article is highly flammable - keep away from heat and/or fire.

donning and doffing	PUTTING ON:
	Place the mask over the mouth and nose on the face and secure it in place with straps or ear loops.  Adjust the flexible nose clip to the shape of the face so that as little air as possible can escape from the upper end.  PULL OUT:  Loosen the straps to remove the face mask from the face. Do not touch the fleece material. Disinfect your hands after disposing of the face mask.

\*slight deviations possible due to standard tolerances

rev-no.: 2025-01 changes and errors excepted