

#### AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

# **Technical Data Sheet**

Article-No.: 05020-W

Description: MED-COMFORT

disposable PE apron light, embossed, in a bag



SIZES	
	Uni
width	75 cm
length	125 cm

PRODUCT DESCRIPTION							
material	☐ Polypro-	✓ Polyethy-	☐ SMS				
	pylene (PP)	lene (PE)					
colour	☐ blue [B]	green [G]	orange [O]	red [R]	white [W]	yellow [Y]	☐ light-blue
							[LB]
	☐ dark-blue	☐ dark-green					
	[DB]	[DG]					
closure	✓ to tie	☐ zipper	☐ push	velcro			
			bottons				
characteristics	☐ tricot-cuff	elastic wrist-	✓ latex-free	☐ hood	☐ stand-up	☐ partially	sewed PE-
		cuff			collar	laminated	oversleeves
		<u> </u>					
	shortenable	☐ blocked	☐ dust-cover		☐ short sleeves	□ with	☐ single folded
	by perforation					cardboard	
	perioration						
REGULATORY	AFFAIRS						
PPE-Regulation	✓ category I	ategory II	☐ category III	no PPE-			
(EU) 2016/425				article			
. , .							CE
MD-Regulation	☐ class I	class II	☐ class III	sterile	☐ measuring	✓ no medical	
(EU) 2017/745					fundtion	device	
Food Contact		□ aqueous		☑ alcoholic	☑ dry foods	not	
(EG) 1935/2004		foods		foods		approved	77
						for food-	, ,
						contact	
STANDARDS							
tested in	☐ DIN EN	☐ ISO 16603	☐ ISO 16604	☐ ISO 22610	☐ ISO 22611	☐ ISO 22612	
accordance to	14126	130 10003	130 10004	130 22010	130 22011	130 22012	
	1.120						
	☐ EN ISO	☐ EN 13034					
	13688						

QMFORM\_67.001 1/3 Herausgabe: 15.09.2021



#### AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

# **Technical Data Sheet**

Article-No.: 05020-W

Description: MED-COMFORT

disposable PE apron light, embossed, in a bag

LOGISTIC DATA SUBPACKING			
general information			
material		bag	
pieces per subpacking	100		
GTIN subpacking		4044941020013	
measures & sizes			
length	480 mm		
width		320 mm	
height -			
weights			
size	net weight	brut weight	
75x125 cm	760 g	765 g	

LOGISTIC DATA PALETTE			
general information			
kind of palette euro-palette			
measures & sizes			
cartons per layer 6			
layer per palette	7		
palette height		151,5	
weights			
size	net weight	brut weight	
uni	342 kg	367 kg	



LOGISTIC DATA OUTER PACKING			
general information			
material carton		carton	
subpackings per outer packing		10	
GTIN outer packing		4044941020211	
measures & sizes			
length		390 mm	
width		300 mm	
height		195 mm	
weights			
size	net weight	brut weight	
75x125 cm	7.650 g	8.150 g	



#### AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

### **Technical Data Sheet**

Article-No.: 05020-W

Description: **MED-COMFORT** 

disposable PE apron light, embossed, in a bag



#### **WARNINGS AND SAFETY INFORMATION**

storage / expiry	Store product in original packaging in a cool and dry place without additional weight, protect from direct sunlight. Do
date	not store near ozone sources (laser printers, copiers). The actual expiry time in use cannot be specified in general terms, as it depends on the general conditions of use. An individual risk assessment must be carried out in each case. The expiry date - valid for proper storage - is stated on the packaging.

use and control	Only ever use the product for its intended use. A check/risk assessment must be carried out to ensure that the product is suitable for the intended use.  Ensure that the size matches that of the user. Do not make any changes to the product.  Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.).  This product is intended for single use only

The used product must be disposed of in accordance with the disposal regulations of the local waste disposal Unused products can be disposed of with household waste.	
warnings	Prolonged wearing can lead to heat stress.  This product is flammable - keep away from fire.

	The product of naminable Reep and y nomine.

#### donning and doffing WAY OF PUTTING ON:

### instructions

- Make sure that the size matches that of the user. Do not make any changes to the product.
- Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.).
- Open the gown and put it on. Make sure that the garment does not tear. Close the gown by wrapping the waistband around your body and fastening it with a bow.

REMOVAL METHOD:

- Undo the waist and collar ties and remove the gown.

rev-no.: 4

21. 10.2024 date:

changes and errors excepted

3/3 Herausgabe: 15.09.2021