

Technical Data Sheet

Article-No.: **02101**
Description: **MED-COMFORT**
facemask
3-ply, to tie



SIZES	
length	width
175 mm	95 mm

PRODUCT DESCRIPTION							
material	☒ non-woven PP	☐ paper					
colour	☒ blue [B]	☐ berry [Berry]	☐ bordeaux [BR]	☐ cedro [C]	☒ green [G]	☐ lagoon [Lagoon]	☐ purple [L]
	☐ lime [Lime]	☐ Mix (blue,yellow,	☐ orange [O]	☐ pink [P]	☒ white [W]	☐ yellow [Y]	☐ platinum / grey [PT]
	☐ blue patterned [DB]	☐ wine-red patterned [WR]	☐ black [SCH]	☐ black-red-gold	☐ clean ocean / ocean green [CO]		
characteristics	☒ type IIR	☒ 3 layers	☒ latex-free	☐ with valve	☐ anti-fog	☒ to tie	☐ with activated carbon filter
characteristics	bacterial filtration capacity ≥ 98%, splash resistance ≥ 16 kPa						

REGULATORY AFFAIRS							
PPE-Regulation (EU) 2016/425	☐ category I	☐ category II	☐ category III	☒ no PPE-article			CE
MD-Regulation (EU) 2017/745	☒ class I	☐ class II	☐ class III	☐ sterile	☐ measuring fundtion	☐ no medical device	

STANDARDS							
tested in accordance to	☒ EN 14683:2019	☐ EN 149:2001					

Technical Data Sheet

Article-No.: **02101**
Description: **MED-COMFORT**
facemask
3-ply, to tie



LOGISTIC DATA SUBPACKING	
general information	
material	carton
pieces per subpacking	50
GTIN subpacking blue [B]	4044941030876
GTIN subpacking berry [B]	-
GTIN subpacking bordeaux [BR]	-
GTIN subpacking cedro [C]	-
GTIN subpacking green [G]	4044941030883
GTIN subpacking lagoon [Lagoon]	-
GTIN subpacking purple [L]	-
GTIN subpacking lime [L]	-
GTIN subpacking mix [Mix]	-
GTIN subpacking orange [O]	-
GTIN subpacking pink [P]	-
GTIN subpacking white [W]	4044941030869
GTIN subpacking yellow [Y]	-
GTIN subpacking platinum [PT]	-
GTIN subpacking blue patterned [DB]	-
GTIN subpacking wine-red patterned [WR]	-
GTIN subpacking black [SCH]	-
GTIN subpacking black-red-gold	-
GTIN subpacking clean ocean [CO]	-
measures & sizes	
length	180 mm
width	110 mm
height	115 mm
weights	
size	brut weight*
uni	250 g

LOGISTIC DATA PALETTE	
general information	
kind of palette	euro-palette
measures & sizes	
cartons per layer	8
layer per palette	7
palette height	190
weights	
size	brut weight*
uni	193 kg

LOGISTIC DATA OUTER PACKING	
general information	
material	carton
subpackings per outer packing	10
GTIN outer packing blue [B]	4044941030937
GTIN outer packing berry [B]	-
GTIN outer packing bordeaux [BR]	-
GTIN outer packing cedro [C]	-
GTIN outer packing green [G]	4044941030944
GTIN outer packing lagoon [Lagoon]	-
GTIN outer packing purple [L]	-
GTIN outer packing lime [L]	-
GTIN outer packing mix [Mix]	-
GTIN outer packing orange [O]	-
GTIN outer packing pink [P]	-
GTIN outer packing white [W]	4044941030920
GTIN outer packing yellow [Y]	-
GTIN outer packing platinum [PT]	-
GTIN outer packing blue patterned [DB]	-
GTIN outer packing wine-red patterned [WR]	-
GTIN outer packing black [SCH]	-
GTIN outer packing black-red-gold	-
GTIN outer packing clean ocean [CO]	-
measures & sizes	
length	570 mm
width	190 mm
height	250 mm
weights	
size	brut weight*
uni	3.000 g

Technical Data Sheet

Article-No.: **02101**
Description: **MED-COMFORT**
facemask
3-ply, to tie



WARNINGS AND SAFETY INFORMATION

storage / expiry date	Store face mask in original packaging in a cool and dry place without additional weight, protect from direct sunlight. Do not store near ozone sources (laser printers, copiers). The actual expiry time in use cannot be specified in general terms, as it depends on the general conditions of use. An individual risk assessment must be carried out in each case. The expiry date - valid for proper storage - is stated on the packaging.
use and control	Always use face masks only for their intended use. A check/risk assessment must be carried out to ensure that the face mask is suitable for the intended use. Change the mask when it is soaked through.
disposal	Used face masks must be disposed of in accordance with the disposal regulations of the local waste disposal company. Unused face masks can be disposed of with household waste.
warnings	Change the mask when it is soaked through. Article is highly flammable - keep away from heat and/or fire.
donning and doffing instructions	PUTTING ON: Place the mask over the mouth and nose on the face and secure it in place with straps or ear loops. Adjust the flexible nose clip to the shape of the face so that as little air as possible can escape from the upper end. PULL OUT: Loosen the straps to remove the face mask from the face. Do not touch the fleece material. Disinfect your hands after disposing of the face mask.

*slight deviations possible due to standard tolerances

rev-no.: 2026-01
date: 09.02.2026

changes and errors excepted