

AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: 05561-Y

Description: **MED-COMFORT**

> non-woven gown with PE-coating



SIZES	
	L/XL
width	138 cm
length	138 cm

EN

PRODUCT DES	RODUCT DESCRIPTION													
material	2	Polypro- pylene (PP)	7	Polyethy- lene (PE)		SMS								
colour		blue [B]		green [G]		orange [O]		red [R]		white [W]	8	yellow [Y]		light-blue [LB]
		dark-blue [DB]		dark-green [DG]										
closure	V	to tie		zipper		push bottons		velcro						
characteristics	v	tricot-cuff		elastic wrist- cuff	7	latex-free		hood		stand-up collar	Z.	partially laminated		sewed PE- oversleeves
		shortenable by perforation		blocked		dust-cover		embossed		short sleeves		with cardboard	V	single folded

REGULATORY AFFAIRS													
PPE-Regulation (EU) 2016/425		category I		category II		category III	>	no PPE- article				((
MD-Regulation (EU) 2017/745	v	class I		class II		class III		sterile	measuring fundtion		no medical device		
Food Contact (EG) 1935/2004		acidic foods		aqueous foods		fatty foods		alcoholic foods	dry foods	7	not approved for food-contact		

STANDARDS												
Material tested	☑ DIN	EN 🗵	ISO 16603	V	ISO 16604	V	ISO 22610	7	ISO 22611	v	ISO 22612	
on the basis of	141	26										
	□ EN I 136	ı	EN 13034									

1/3 Herausgabe: 15.09.2021



AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: **05561-Y**

Description: MED-COMFORT

non-woven gown with PE-coating

LOGISTIC DATA SUBPACKING							
general information							
material		bag					
pieces per subpacking	5	10					
GTIN subpacking blue	[B]	-					
GTIN subpacking gree	en [G]	-					
GTIN subpacking orar	nge [O]	-					
GTIN subpacking red	[R]	-					
GTIN subpacking whit	-						
GTIN subpacking yello	ow [Y]	4044941022451					
GTIN subpacking light	-blue [LB]	-					
GTIN subpacking dark	c-blue [DB]	-					
GTIN subpacking dark	-green [DG]	-					
measures & sizes							
length		430 mm					
width	330 mm						
height	-						
weights							
size	net weight	brut weight					
L/XL	650 g	655 g					

LOGISTIC DATA PALETTE							
general information							
kind of palette		euro-palette					
measures & sizes							
cartons per layer		4					
layer per palette		6					
palette height		201					
weights	weights						
size	net weight	brut weight					
L/XL	169 kg	194 kg					

EN



LOGISTIC DATA OUTER PACKING							
general information							
material		carton					
subpackings per outer	packing	10					
GTIN outer packing blu	ue [B]	-					
GTIN outer packing gr	een [G]	-					
GTIN outer packing or	ange [O]	-					
GTIN outer packing re	d [R]	-					
GTIN outer packing wh	nite [W]	-					
GTIN outer packing ye	llow [Y]	4044941022468					
GTIN outer packing lig	ht-blue [LB]	-					
GTIN outer packing da	rk-blue [DB]	-					
GTIN outer packing da	rk-green [DG]	-					
measures & sizes							
length		500 mm					
width		280 mm					
height		310 mm					
weights							
size	net weight	brut weight					
L/XL	6.550 g	7.050 g					

QMFORM_67.001 2/3 Herausgabe: 15.09.2021



AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: 05561-Y

Description: **MED-COMFORT**

> non-woven gown with PE-coating



WARNINGS AND SAFETY INFORMATION

storage / expiry	Store product in original packaging in a cool and dry place without additional weight, protect from direct sunlight. Do
date	not store near ozone sources (laser printers, copiers). The actual expiry time in use cannot be specified in general terms,
	as it depends on the general conditions of use. An individual risk assessment must be carried out in each case. The
	expiry date - valid for proper storage - is stated on the packaging.

use and control	Only ever use the product for its intended use. A check/risk assessment must be carried out to ensure that the product	l
	is suitable for the intended use.	1
	Ensure that the size matches that of the user. Do not make any changes to the product.	l
	Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.).	l
	This product is intended for single use only	l
		l
		l

disposal	The used product must be disposed of in accordance with the disposal regulations of the local waste disposal company.
	Unused products can be disposed of with household waste.

warnings	Prolonged wearing can lead to heat stress.	
	This product is flammable - keep away from fire.	

donning and doffing WAY OF PUTTING ON:

instructions

- Make sure that the size matches that of the user. Do not make any changes to the product.
- Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.).

3/3

- Open the gown and put it on. Make sure that the garment does not tear. Close the gown by wrapping the waistband around your body and fastening it with a bow.

REMOVAL METHOD:

- Undo the waist and collar ties and remove the gown.

rev-no.:

date: 04. 12.2024 changes and errors excepted

QMFORM_67.001 Herausgabe: 15.09.2021