

AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: 05561-Y

Description: **MED-COMFORT**

> non-woven gown with PE-coating



| SIZES | |
|--------|--------|
| | L/XL |
| width | 138 cm |
| length | 138 cm |

EN

| PRODUCT DES | CRI | IPTION | | | | | | | | | | |
|--------------------|-----|----------------------------------|----|------------------------|---|-----------------|----------|--------------------|----|------------------------|---|--------------------------|
| material | Ø | Polypro- pylene (PP) | Į. | Polyethy- lene (PE) | | SMS | | | | | | |
| colour | | blue [B] | | green [G] | | orange [O] | red [R] | white [W] | Z | yellow [Y] | | light-blue [LB] |
| | | dark-blue [DB] | | dark-green [DG] | | | | | | | | |
| closure | | to tie | | zipper | | push bottons | velcro | | | | | |
| characteristics | v | tricot-cuff | | elastic wrist- cuff | 7 | latex-free | hood | stand-up collar | Z. | partially laminated | | sewed PE- oversleeves |
| | | shortenable by perforation | | blocked | | dust-cover | embossed | short sleeves | | with cardboard | V | single folded |

| REGULATORY | ΑFI | FAIRS | | | | | | | | |
|---------------------------------|-----|--------------|------------------|--------------|---|--------------------|--------------------|---|-------------------------------------|----|
| PPE-Regulation (EU) 2016/425 | | category I | category II | category III | 7 | no PPE- article | | | | (6 |
| MD-Regulation (EU) 2017/745 | v | class I | class II | class III | | sterile | measuring fundtion | | no medical device | |
| Food Contact (EG) 1935/2004 | | acidic foods | aqueous foods | fatty foods | | alcoholic foods | dry foods | V | not approved for food-contact | |

| STANDARDS | | | | | | | |
|---------------|-------------------|-------------|-------------|-------------|-------------|-------------|--|
| tested in | □ DIN EN | □ ISO 16603 | □ ISO 16604 | □ ISO 22610 | □ ISO 22611 | □ ISO 22612 | |
| accordance to | 14126 | | | | | | |
| | □ EN ISO 13688 | □ EN 13034 | | | | | |

1/3 Herausgabe: 15.09.2021



AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany


Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: 05561-Y

Description: MED-COMFORT

non-woven gown with PE-coating



| LOGISTIC DATA PA | LOGISTIC DATA PALETTE | | | | | | |
|---------------------|-----------------------|-------------|--|--|--|--|--|
| general information | | | | | | | |
| kind of palette | euro-palette | | | | | | |
| measures & sizes | measures & sizes | | | | | | |
| cartons per layer | 4 | | | | | | |
| layer per palette | | 6 | | | | | |
| palette height | 201 | | | | | | |
| weights | weights | | | | | | |
| size | net weight | brut weight | | | | | |
| L/XL | 169 g | 194 g | | | | | |



| LOGISTIC DATA OU | TER PACKING | |
|------------------------|---------------|-------------|
| general information | | |
| material | carton | |
| subpackings per outer | packing | 10 |
| GTIN outer packing blu | ue [B] | - |
| GTIN outer packing gr | een [G] | - |
| GTIN outer packing or | ange [O] | - |
| GTIN outer packing re- | d [R] | - |
| GTIN outer packing wh | - | |
| GTIN outer packing ye | 4044941022468 | |
| GTIN outer packing lig | ht-blue [LB] | - |
| GTIN outer packing da | rk-blue [DB] | - |
| GTIN outer packing da | rk-green [DG] | - |
| measures & sizes | | |
| length | | 500 mm |
| width | | 280 mm |
| height | 310 mm | |
| weights | | |
| size | net weight | brut weight |
| L/XL | 6.550 g | 7.050 g |



AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: 05561-Y

Description: **MED-COMFORT**

> non-woven gown with PE-coating



WARNINGS AND SAFETY INFORMATION

| storage / expiry | Store product in original packaging in a cool and dry place without additional weight, protect from direct sunlight. Do |
|------------------|---|
| date | not store near ozone sources (laser printers, copiers). The actual expiry time in use cannot be specified in general terms, |
| | as it depends on the general conditions of use. An individual risk assessment must be carried out in each case. The |
| | expiry date - valid for proper storage - is stated on the packaging. |
| | |

| use and control | Only ever use the product for its intended use. A check/risk assessment must be carried out to ensure that the product | |
|-----------------|--|--|
| | is suitable for the intended use. | |
| | Ensure that the size matches that of the user. Do not make any changes to the product. | |
| | Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.). | |
| | This product is intended for single use only | |
| | | |
| | | |

| disposal | The used product must be disposed of in accordance with the disposal regulations of the local waste disposal company. |
|----------|---|
| | Unused products can be disposed of with household waste. |

| warnings | Prolonged wearing can lead to heat stress. | |
|----------|--|--|
| | This product is flammable - keep away from fire. | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |

donning and doffing WAY OF PUTTING ON:

instructions

- Make sure that the size matches that of the user. Do not make any changes to the product.
- Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.).

3/3

- Open the gown and put it on. Make sure that the garment does not tear. Close the gown by wrapping the waistband around your body and fastening it with a bow.

REMOVAL METHOD:

- Undo the waist and collar ties and remove the gown.

rev-no.:

date: 04. 12.2024 changes and errors excepted

QMFORM_67.001 Herausgabe: 15.09.2021