

Technical Data Sheet

Article-No.: **09820**
 Description: **MED-COMFORT**
tourniquet
Individually packed



DIMENSIONS				
length	width	height	diameter	capacity
450 mm	25 mm	-	-	-

PRODUCT DESCRIPTION							
material	☐ Polypropylene (PP)	☐ Polyethylene (PE)	☐ PP non woven	☐ PP non woven with PE coating	☐ paper	☒ Clip: ABS Band: latex-free elastic strap	
colour	☒ blue [B]	☐ cedro [C]	☐ dark blue [DB]	☒ green [G]	☐ green opac [GO]	☐ purple [L]	☐ lightblue [LB]
	☐ Lime [LIME]	☐ Mix [MIX] black, cedro, purple, pink, lightblue	☐ orange [O]	☐ pink [P]	☒ red [R]	☐ rose [ROSE]	☐ black [SCH]
	☐ violett [V]	☒ white [W]	☐ wine red [WR]	☒ yellow [Y]	☐ silver	☐ natural	☐ bordeaux
characteristics	☐ single use	☐ single bladed	☐ 5 ply	☐ lint free	☐ 6 ply	☒ elastic	
ingredients	☐ 50% Urea (CAS 57-13-6)	☐ 50% water	☐ 70 % isopropyl alcohol				
remark							

REGULATORY AFFAIRS							
PPE-Regulation (EU) 2016/425	☐ Category I	☐ Category II	☐ Category III	☒ no PPE-article			CE
#BEZUG!	☒ Class I	☐ Class II	☐ Class III	☐ sterile	☐ measuring function	☐ no medical device	
Food Contact (EG) 1935/2004	☐ acidic foods	☐ aqueous foods	☐ fatty foods	☐ alcoholic foods	☐ dry foods	☒ not approved for food-contact	

Technical Data Sheet

Article-No.: **09820**
Description: **MED-COMFORT**
tourniquet
Individually packed



LOGISTIC DATA SUBPACKING	
general information	
material	carton/PVC
pieces per subpacking	1
GTIN subpacking blue [B]	4044941016245
GTIN subpacking cedro [C]	-
GTIN subpacking dark blue [DB]	-
GTIN subpacking green [G]	4044941010359
GTIN subpacking green opac [GO]	-
GTIN subpacking purple [L]	-
GTIN subpacking lightblue [LB]	-
GTIN subpacking Lime [LIME]	-
GTIN subpacking Mix [MIX]	-
GTIN subpacking orange [O]	4044941727196
GTIN subpacking pink [P]	-
GTIN subpacking rot [R]	4044941010335
GTIN subpacking rose [ROSE]	-
GTIN subpacking black [SCH]	-
GTIN subpacking violet [V]	-
GTIN subpacking white [W]	-
GTIN subpacking wine red [WR]	-
GTIN subpacking yellow [Y]	4044941016269
GTIN subpacking silver	-
GTIN subpacking natur	-
GTIN subpacking	-
measures & size	
length	42 mm
width	40 mm
height	85 mm
weights	
size	gross weight*
uni	36 g

LOGISTIC DATA OUTER PACKING	
general information	
material	carton
subpackings per outer packing	400
GTIN oute packing blue [B]	4044941016252
GTIN oute packing cedro [C]	-
GTIN oute packing dark blue [DB]	-
GTIN oute packing green [G]	4044941010366
GTIN oute packing green opac [GO]	-
GTIN oute packing purple [L]	-
GTIN oute packing lightblue [LB]	-
GTIN oute packing Lime [LIME]	-
GTIN oute packing Mix [MIX]	-
GTIN oute packing orange [O]	4044941727202
GTIN oute packing pink [P]	-
GTIN oute packing rot [R]	4044941010342
GTIN oute packing rose [ROSE]	-
GTIN oute packing black [SCH]	-
GTIN oute packing violet [V]	-
GTIN oute packing white [W]	-
GTIN oute packing wine red [WR]	-
GTIN oute packing yellow [Y]	4044941016276
GTIN oute packing silver	-
GTIN oute packing natur	-
GTIN oute packing	-
measures & size	
length	453 mm
width	380 mm
height	423 mm
weights	
size	gross weight*
uni	14.900 g

LOGISTIC DATA PALETTE	
general information	
kind of palette	euro-palette
measures & size	
cartons per layer	4
layer per palette	4
height of the palette	184 cm
weights	
size	brut weight*
uni	263 kg

Technical Data Sheet

Article-No.: **09820**
Description: **MED-COMFORT**
tourniquet
Individually packed



WARNINGS AND SAFETY INFORMATION

storage	Store items in original packaging in a cool and dry place without additional weight, protect from direct sunlight.
use and control	1. Open the package and take out the tourniquet, expand it with the face of the fastener upward, then apply the stretch strap around appropriate area or cubitus area. 2. Insert the fastener into the insertion part and fasten it. 3. Pull the stretch strap till it reaches appropriate tightness. 4. While releasing the tourniquet, first press the face of the fastener to release the strap. Then press the button on the face, the insertion part will leave the fastener automatically.
cleaning instructions	With low contamination, hand washing is possible at max. 60 °C with household washing detergent. Disinfection with liquid disinfectant is recommended. Not suitable for physical disinfection procedures.
disposal	The product must be disposed of in accordance with the disposal regulations of the local waste disposal company. Unused products can be disposed of with household waste. Do not dispose of in the toilet.
warnings / safety information	Use only by trained medical personnel. Improper use can lead to injuries or complications.

rev-no.: 2025-01
date: 09. 09.2025

*slight deviations possible due to standard tolerances

changes and errors excepted