

AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: **02105-B**

SIZES length

EN

Description: MED-COMFORT

width

facemask 4-ply, to tie with fog-stop



175 mm	95 mm						
PRODUCT DES	SCRIPTION						
material	☑ non-woven	☐ paper					
colour	☑ blue [B]	☐ berry [Berry]	□ bordeaux [BR]	☐ cedro [C]	☐ green [G]	☐ lagoon [Lagoon]	☐ purple [L]
	☐ lime [Lime]	☐ Mix (blue,yellow,	☐ orange [O]	☐ pink [P]	☐ white [W]	☐ yellow [Y]	☐ platinum / grey [PT]
	blue patterned [DB]	wine-red patterned [WR]	□ black [SCH]	black-red- gold	clean ocean / ocean green [CO]		
characteristics	✓ type IIR	✓ 4 layers	☑ latex-free	☐ with valve	☑ anti-fog	☑ to tie	☐ withactivatedcarbon filter
characteristics	bacterial filtratio	n capacity ≥ 98%, s	plash resistance :	≥ 16 kPa			
REGULATORY	AFFAIRS						
PPE-Regulation (EU) 2016/425	☐ category I	☐ category II	☐ category III	ono PPE- article			CE
MD-Regulation (EU) 2017/745	✓ class I	☐ class II	☐ class III	sterile	measuring fundtion	no medical device	
STANDARDS							
tested in accordance to	✓ EN 14683:2019	☐ EN 149:2001					

QMFORM_67.001 1/3 Herausgabe: 15.09.2021



AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany


Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: 02105-B

Description: MED-COMFORT

facemask 4-ply, to tie with fog-stop



LOGISTIC DATA PALETTE	
general information	
kind of palette	euro-palette
measures & sizes	
cartons per layer	7
layer per palette	6
palette height	183
weights	
size	brut weight*
uni	166 kg





LOGISTIC DATA OUTER PACKING	
general information	
material	carton
subpackings per outer packing	10
GTIN outer packing blue [B]	4044941730592
GTIN outer packing berry [B]	-
GTIN outer packing bordeaux [BR]	-
GTIN outer packing cedro [C]	-
GTIN outer packing green [G]	-
GTIN outer packing lagoon [Lagoon]	-
GTIN outer packing purple [L]	-
GTIN outer packing lime [L]	-
GTIN outer packing mix [Mix]	-
GTIN outer packing orange [O]	-
GTIN outer packing pink [P]	-
GTIN outer packing white [W]	-
GTIN outer packing yellow [Y]	-
GTIN outer packing platinum [PT]	-
GTIN outer packing blue patterned [DB]	-
GTIN outer packing wine-red patterned [WF	-
GTIN outer packing black [SCH]	-
GTIN outer packing black-red-gold	-
GTIN outer packing clean ocean [CO]	-
measures & sizes	
length	570 mm
width	205 mm
height	280 mm
weights	
size	brut weight*
uni	3.350 g



AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: 02105-B

Description: **MED-COMFORT**

> facemask 4-ply, to tie with fog-stop





WARNINGS AND SAFETY INFORMATION

storage / expiry	Store face mask in original packaging in a cool and dry place without additional weight, protect from direct sunlight. Do
date	not store near ozone sources (laser printers, copiers). The actual expiry time in use cannot be specified in general terms, as it depends on the general conditions of use. An individual risk assessment must be carried out in each case. The expiry date - valid for proper storage - is stated on the packaging.

use and control	Always use face masks only for their intended use. A check/risk assessment must be carried out to ensure that the face
	mask is suitable for the intended use.
	Change the mask when it is soaked through.

disposal	Used face masks must be disposed of in accordance with the disposal regulations of the local waste disposal company.
	Unused face masks can be disposed of with household waste.

warnings	Change the mask when it is soaked through. Article is highly flammable - keep away from heat and/or fire.

donning and doffing	PUTTING ON:
	Place the mask over the mouth and nose on the face and secure it in place with straps or ear loops. Adjust the flexible nose clip to the shape of the face so that as little air as possible can escape from the upper end. PULL OUT: Loosen the straps to remove the face mask from the face. Do not touch the fleece material. Disinfect your hands after disposing of the face mask.

*slight deviations possible due to standard tolerances

rev-no.: 2025-01 17. 11.2025 changes and errors excepted date:

EN 3/3 Herausgabe: 15.09.2021