

Technical Data Sheet

Article-No.: **05020-BS-Block-140**
Description: **MED-COMFORT**
disposable PE apron
light, smooth, blocked, in a bag



SIZES	
	Uni
width	75 cm
length	140 cm

PRODUCT DESCRIPTION							
material	<input type="checkbox"/> Polypropylene (PP)	<input checked="" type="checkbox"/> Polyethylene (PE)	<input type="checkbox"/> SMS				
colour	<input checked="" type="checkbox"/> blue [B]	<input type="checkbox"/> green [G]	<input type="checkbox"/> orange [O]	<input type="checkbox"/> red [R]	<input type="checkbox"/> white [W]	<input type="checkbox"/> yellow [Y]	<input type="checkbox"/> light-blue [LB]
	<input type="checkbox"/> dark-blue [DB]	<input type="checkbox"/> dark-green [DG]					
closure	<input checked="" type="checkbox"/> to tie	<input type="checkbox"/> zipper	<input type="checkbox"/> push buttons	<input type="checkbox"/> velcro			
characteristics	<input type="checkbox"/> tricot-cuff	<input type="checkbox"/> elastic wrist-cuff	<input checked="" type="checkbox"/> latex-free	<input type="checkbox"/> hood	<input type="checkbox"/> stand-up collar	<input type="checkbox"/> partially laminated	<input type="checkbox"/> sewed PE-oversleeves
	<input type="checkbox"/> shortenable by perforation	<input checked="" type="checkbox"/> blocked	<input type="checkbox"/> dust-cover	<input type="checkbox"/> embossed	<input type="checkbox"/> short sleeves	<input type="checkbox"/> with cardboard	<input type="checkbox"/> single folded

REGULATORY AFFAIRS							
PPE-Regulation (EU) 2016/425	<input checked="" type="checkbox"/> category I	<input type="checkbox"/> category II	<input type="checkbox"/> category III	<input type="checkbox"/> no PPE-article			CE
MD-Regulation (EU) 2017/745	<input type="checkbox"/> class I	<input type="checkbox"/> class II	<input type="checkbox"/> class III	<input type="checkbox"/> sterile	<input type="checkbox"/> measuring function	<input checked="" type="checkbox"/> no medical device	
Food Contact (EG) 1935/2004	<input checked="" type="checkbox"/> acidic foods	<input checked="" type="checkbox"/> aqueous foods	<input checked="" type="checkbox"/> fatty foods	<input checked="" type="checkbox"/> alcoholic foods	<input checked="" type="checkbox"/> dry foods	<input type="checkbox"/> not approved for food-contact	☞

STANDARDS							
tested in accordance to	<input type="checkbox"/> DIN EN 14126	<input type="checkbox"/> ISO 16603	<input type="checkbox"/> ISO 16604	<input type="checkbox"/> ISO 22610	<input type="checkbox"/> ISO 22611	<input type="checkbox"/> ISO 22612	
	<input type="checkbox"/> EN ISO 13688	<input type="checkbox"/> EN 13034					

Technical Data Sheet

Article-No.: **05020-BS-Block-140**
 Description: **MED-COMFORT**
disposable PE apron
light, smooth, blocked, in a bag



LOGISTIC DATA SUBPACKING		
general information		
material	bag	
pieces per subpacking	50	
GTIN subpacking	4044941022178	
measures & sizes		
length	480 mm	
width	320 mm	
height	-	
weights		
size	net weight	brut weight
75x140 cm	660 g	665 g

LOGISTIC DATA PALETTE		
general information		
kind of palette	euro-palette	
measures & sizes		
cartons per layer	6	
layer per palette	10	
palette height	150	
weights		
size	net weight	brut weight
75x140 cm	429 g	454 g

LOGISTIC DATA OUTER PACKING		
general information		
material	carton	
subpackings per outer packing	10	
GTIN outer packing	4044941022185	
measures & sizes		
length	390 mm	
width	300 mm	
height	135 mm	
weights		
size	net weight	brut weight
75x140 cm	6.650 g	7.150 g

Technical Data Sheet

Article-No.: **05020-BS-Block-140**
Description: **MED-COMFORT**
disposable PE apron
light, smooth, blocked, in a bag



WARNINGS AND SAFETY INFORMATION

storage / expiry date	Store product in original packaging in a cool and dry place without additional weight, protect from direct sunlight. Do not store near ozone sources (laser printers, copiers). The actual expiry time in use cannot be specified in general terms, as it depends on the general conditions of use. An individual risk assessment must be carried out in each case. The expiry date - valid for proper storage - is stated on the packaging.
use and control	Only ever use the product for its intended use. A check/risk assessment must be carried out to ensure that the product is suitable for the intended use. Ensure that the size matches that of the user. Do not make any changes to the product. Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.). This product is intended for single use only
disposal	The used product must be disposed of in accordance with the disposal regulations of the local waste disposal company. Unused products can be disposed of with household waste.
warnings	Prolonged wearing can lead to heat stress. This product is flammable - keep away from fire.
donning and doffing instructions	WAY OF PUTTING ON: - Make sure that the size matches that of the user. Do not make any changes to the product. - Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.). - Open the gown and put it on. Make sure that the garment does not tear. Close the gown by wrapping the waistband around your body and fastening it with a bow. REMOVAL METHOD: - Undo the waist and collar ties and remove the gown.

rev-no.: 4
date: 21.10.2024

changes and errors excepted