

Technical Data Sheet

Article-No.: **09097**

Description: **MED-COMFORT**

bed protection underlay

reusable



DIMENSIONS	
length	750 cm
width	900 cm
height	-

PRODUCT DESCRIPTION							
material	☐ polypropylene (PP)	☒ polyethylene (PE)	☐ PP fleece with PE coating	☐ cellulose with PE film	☒ 1st layer: Polyester, white 2nd layer: Absorbent core made of viscose 3rd layer: Liquid barrier made of polyurethane 4th layer: Polyester, blue		
colour	☒ white	☒ blue	☐ green	☐ yellow	☐ red	☐ black	
characteristics	☐ embossed	☐ 2 layers	☐ 8 threads	☒ washable up to 90°C	☒ reusable		
suction volume/absorption	☐ ca. 200 ml	☐ ca. 525 ml	☐ ca. 780 ml	☐ 382 ml ± 10%	☐ 915 ml ± 10%	☐ 1.200 ml ± 10%	☒ 1.300 ml ± 10%
perforation	☐ 38 cm	☐ 50 cm					

REGULATORY AFFAIRS							
PPE-Regulation (EU) 2016/425	☐ Category I	☐ Category II	☐ Category III	☒ no PPE-article			CE
MD-Regulation (EU) 2017/745	☒ Class I	☐ Class II	☐ Class III	☐ sterile	☐ measuring function	☐ no medical device	
Food Contact (EG) 1935/2004	☐ acidic foods	☐ aqueous foods	☐ fatty foods	☐ alcoholic foods	☐ dry foods	☒ not approved for food-contact	

LOGISTIC DATA SUBPACKING	
general information	
material	bag
pieces per subpacking	1
EAN subpacking	4044941756394
measures & size	
length	350 mm
width	320 mm
height	-
weights	
size	brut weight*
uni	335 g
PZN/PPN	
PZN	20048894
PPN	112004889436

LOGISTIC DATA OUTER PACKING	
general information	
material	carton
subpackings per outer packing	20
EAN outer packing blue	4044941756400
measures & size	
length	520 mm
width	310 mm
height	290 mm
weights	
size	brut weight*
uni	7.200 g

LOGISTIC DATA PALETTE	
general information	
kind of palette	euro-palette
measures & size	
cartons per layer	4
layer per palette	6
height of the palette	189 cm
weights	
size	brut weight*
uni	198 kg

Technical Data Sheet

Article-No.: **09097**
Description: **MED-COMFORT**
bed protection underlay
reusable



WARNINGS AND SAFETY INFORMATION

storage	Store items in original packaging in a cool and dry place without additional weight, protect from direct sunlight.
use and control	Place the blue side face down on the couch/stretcher and the white side facing up to ensure optimal fluid absorption. Check the product for proper condition before use and replace it if necessary.
cleaning instructions	Machine washable up to 90 degrees, suitable for tumble drying, do not iron, do not bleach, do not dry clean.
disposal	The used product must be disposed of in accordance with the disposal regulations of the local waste disposal company. Unused products can be disposed of with household waste.
warnings / safety information	Please observe the care instructions. For external use only and only for the intended purpose. Do not use if the product is damaged.

*slight deviations possible due to standard tolerances

rev-no.: 2025-03
date: 13. 11.2025

changes and errors excepted