

## Technical Data Sheet

Article-No.: **09280-ODK**

Description: **Holder for angular urine bottles**  
Without cap and chain



SIZES		
length	width	height
310 mm	135 mm	-

PRODUCT DESCRIPTION							
material	<input type="checkbox"/> stainless steel	<input type="checkbox"/> acrylic glass	<input checked="" type="checkbox"/> PVC coated wire	<input type="checkbox"/> ABS plastic	<input type="checkbox"/> casing: aluminium	<input type="checkbox"/> pump: plastic	<input type="checkbox"/> Polypropylene (PP)
colour	<input type="checkbox"/> silver	<input type="checkbox"/> transparent	<input checked="" type="checkbox"/> white	<input type="checkbox"/> blue			
characteristics	<input type="checkbox"/> with movable hooks	<input type="checkbox"/> with fixing material	<input type="checkbox"/> adjustable stop for depth adjustment	<input checked="" type="checkbox"/> to attach			
desription							
suitable for	all square urin bottles						

LOGISTIC DATA SUBPACKING		
general information		
material	bag	
pieces per subpacking	1	
GTIN subpacking	-	
measures & sizes		
length	500 mm	
width	250 mm	
height	-	
weights		
size	net weight	brut weight
uni	185 g	190 g

LOGISTIC DATA OUTER PACKING		
general information		
material	carton	
subpackings per outer packing	25	
GTIN outer packing	-	
measures & sizes		
length	640 mm	
width	375 mm	
height	395 mm	
weights		
size	net weight	brut weight
uni	4.750 g	5.250 g

LOGISTIC DATA PALETTE		
general information		
kind of palette	euro-palette	
measures & sizes		
cartons per layer	3	
layer per palette	5	
palette height	213 cm	
weights		
size	net weight	brut weight
uni	79 kg	104 kg

## Technical Data Sheet

Article-No.: **09280-ODK**

Description: **Holder for angular urine bottles**  
**Without cap and chain**



### WARNINGS AND SAFETY INFORMATION

<b>storage</b>	Store items in original packaging in a cool and dry place without additional weight, protect from direct sunlight.
<b>warnings / safety information</b>	<ul style="list-style-type: none"><li>- The dispenser is designed exclusively for storing square urine bottles.</li><li>- Product can be hung on the bedside</li><li>- Before use, ensure that all parts are securely mounted and intact.</li></ul>
<b>cleaning instructions</b>	Can be cleaned easily. Use suitable cleaning agents and soft cloths to avoid scratches.