


Technical Data Sheet

Article-No.: **09287**

Description: **urine bag holder**



SIZES		
length	width	height
150 mm	148 mm	-

PRODUCT DESCRIPTION							
material	<input type="checkbox"/> stainless steel	<input type="checkbox"/> acrylic glass	<input checked="" type="checkbox"/> PVC coated wire	<input type="checkbox"/> ABS plastic	<input type="checkbox"/> casing: aluminium	<input type="checkbox"/> pump: plastic	<input type="checkbox"/> Polypropylene (PP)
colour	<input type="checkbox"/> silver	<input type="checkbox"/> transparent	<input checked="" type="checkbox"/> white	<input type="checkbox"/> blue			
characteristics	<input type="checkbox"/> with movable hooks	<input type="checkbox"/> with fixing material	<input type="checkbox"/> adjustable stop for depth adjustment		<input checked="" type="checkbox"/> to attach		
description	to attach						
suitable for							

LOGISTIC DATA SUBPACKING		
general information		
material	bag	
pieces per subpacking	10	
GTIN subpacking	4044941010373	
measures & sizes		
length	500 mm	
width	250 mm	
height	-	
weights		
size	net weight	brut weight
uni	435 g	440 g


LOGISTIC DATA OUTER PACKING		
general information		
material	carton	
subpackings per outer packing	15	
GTIN outer packing	4044941010380	
measures & sizes		
length	510 mm	
width	250 mm	
height	340 mm	
weights		
size	net weight	brut weight
uni	6.600 g	7.100 g

LOGISTIC DATA PALETTE		
general information		
kind of palette	euro-palette	
measures & sizes		
cartons per layer	6	
layer per palette	6	
palette height	219 cm	
weights		
size	net weight	brut weight
uni	256 kg	281 kg

Technical Data Sheet

Article-No.: **09287**

Description: **urine bag holder**



WARNINGS AND SAFETY INFORMATION

storage	Store items in original packaging in a cool and dry place without additional weight, protect from direct sunlight.
warnings / safety information	<ul style="list-style-type: none">- The dispenser is intended exclusively for the storage of urine bags.- Product can be hung on the bedside- Before use, ensure that all parts are securely mounted and intact.
cleaning instructions	Can be cleaned easily. Use suitable cleaning agents and soft cloths to avoid scratches.