

Technical Data Sheet

Article-No.: **09053**
 Description: **Stainless steel holder**
 for PE gloves



SIZES		
length	width	height
165 mm	195 mm	18 mm

PRODUCT DESCRIPTION							
material	<input checked="" type="checkbox"/> stainless steel	<input type="checkbox"/> acrylic glass	<input type="checkbox"/> PVC coated wire	<input type="checkbox"/> ABS plastic	<input type="checkbox"/> casing: aluminium	<input type="checkbox"/> pump: plastic	<input type="checkbox"/> polymethylmethacrylat (PMMA)
colour	<input checked="" type="checkbox"/> silver	<input type="checkbox"/> transparent	<input type="checkbox"/> white	<input type="checkbox"/> blue			
characteristics	<input type="checkbox"/> with movable hooks	<input checked="" type="checkbox"/> with fixing material	<input type="checkbox"/> adjustable stop for depth adjustment	<input checked="" type="checkbox"/> for wall mounting			
description							
suitable for	see dispenser matrix						

LOGISTIC DATA SUBPACKING		
general information		
material	bag	
pieces per subpacking	1	
GTIN subpacking	4044941728513	
measures & sizes		
length	235 mm	
width	170 mm	
height	17 mm	
weights		
size	net weight	brut weight
uni	520 g	525 g

LOGISTIC DATA OUTER PACKING		
general information		
material	carton	
subpackings per outer packing	10	
GTIN outer packing	4044941728520	
measures & sizes		
length	580 mm	
width	440 mm	
height	215 mm	
weights		
size	net weight	brut weight
uni	5.250 g	5.750 g

LOGISTIC DATA PALETTE		
general information		
kind of palette	euro-palette	
measures & sizes		
cartons per layer	2	
layer per palette	9	
palette height	209 cm	
weights		
size	net weight	brut weight
uni	104 kg	129 kg

Technical Data Sheet

Article-No.: **09053**

Description: **Stainless steel holder
for PE gloves**



WARNINGS AND SAFETY INFORMATION

storage	Store items in original packaging in a cool and dry place without additional weight, protect from direct sunlight.
warnings / safety information	<ul style="list-style-type: none">- The dispenser is intended exclusively for the storage and removal of glove boxes.- Mount the product firmly and securely on the wall. Check regularly that the bracket is stable.- When installing the dispenser, ensure that the walls provide the necessary stability- Before use, ensure that all parts are securely mounted and intact.
cleaning instructions	Can be cleaned easily. Use suitable cleaning agents and soft cloths to avoid scratches.