AMPri Handelsgesellschaft mbH
Benzstraße 16
21423 Winsen (Luhe) - Germany
Telefon: +49 (0)4171 / 8480-0
Homepage: www.ampri.de
e-mail: info@ampri.de

Technical Data Sheet



Article-No.: **01240**

Description: **MED-COMFORT Blue Ultra 400**
Nitrile examination glove
blue, non sterile, powder free

THICKNESS		
cuff	**palm**	**fingertips**
0,09 mm	0,15 mm	0,20 mm

PRODUCT DESCRIPTION							
material	☐ Latex	☑ Nitrile	☐ Vinyl	☐ Vinyl-Nitrile-mixture	☐ Polyethy-lene (PE)	☐ TPE	☐ cotton
colour	☐ white	☑ blue	☐ black	☐ mint	☐ purple	☐ mix	☐ cobalt-blue
characteristics	☐ prepowdered	☑ powderfree	☐ sterile	☑ non sterile	☑ ambidex-trous	☐ fits hand-specific	☐ biodgra-dable
surface	☑ full textured	☐ finger textured	☐ not textured	☐ embossed	☑ chlorinated inside		

SIZES							
	XS (5-6)	**S (6-7)**	**M (7-8)**	**L (8-9)**	**XL (9-10)**	**XXL (10-11)**	**XXXL (11-12)**
width	-	80 ± 10 mm	95 ± 10 mm	110 ± 10 mm	115 ± 10 mm	125 ± 10 mm	-
length	-	≥ 400 mm	≥ 400 mm	≥ 400 mm	≥ 400 mm	≥ 400 mm	-

REGULATORY AFFAIRS							
PPE-Regulation (EU) 2016/425	☐ Category I	☐ Category II	☑ Category III	☐ no PPE-article			CE
MD-Regulation (EU) 2017/745	☑ Class I	☐ Class II a	☐ Class III	☐ sterile	☐ measuring function	☐ no medical device	
Food Contact (EG) 1935/2004	☑ acidic foods	☑ aqueous foods	☑ fatty foods	☑ alcoholic foods	☑ dry foods	☐ not approved for food-contact	

STANDARDISATION							
EN 388 Mechanical Risks	**abrasion resistance**	**blade cut resistance Coupe-Test**	**tear resistance**	**puncture resistance**	**blade cut resistance TDM-Test**	**impact test**	
Level	not applicable						

EN 374-1 Chemical Risks	**chemical**	**code letter**	**level**	**permeation time**	**degradation**	ISO 374-1/Type B
	Sodium hydroxide 40%	K	6	> 480 min	-5,8 %	
	Hydrogen Peroxide 30%	P	3	> 60 min	17,2 %	
EN 374-4 Degradation	Formaldehyde 37%	T	6	> 480 min	26,1 %	KPT
EN 374-5 microorganism tightness	The glove is tight against microorganisms (viral, bacteria and fungi). Test according to ISO 16604 - method B.					EN ISO 374-5:2016 VIRUS
EN 420 protective gloves	The glove meets the requirements according to EN 420					
EN 455 medical gloves for single use	The glove meets the requirements according to EN 455-1, EN 455-2, EN 455-3, EN 455-4					EN 455
EN 455-1 freedom from holes	The glove has an AQL < 1.5 in regards to the water retention test (sampling inspection in acc. to ISO 2859-1, general Inspection Level I)					AQL 1.5
ISO 18889 Pesticide application	not applicable					

AMPri Handelsgesellschaft mbH
Benzstraße 16
21423 Winsen (Luhe) - Germany
Telefon: +49 (0)4171 / 8480-0
Homepage: www.ampri.de
e-mail: info@ampri.de

Technical Data Sheet



Article-No.: **01240**

Description: **MED-COMFORT Blue Ultra 400**
Nitrile examination glove
blue, non sterile, powder free

LOGISTIC DATA SUBPACKING	
generell information	
material	carton
pieces per subpacking	50
GTIN subpacking size XS	-
GTIN subpacking size S	4044941007250
GTIN subpacking size M	4044941007267
GTIN subpacking size L	4044941007274
GTIN subpacking size XL	4044941007281
GTIN subpacking size XXL	4044941014180
GTIN subpacking size XXXL	-
PZN subpacking size XS	-
PZN subpacking size S	15318400
PZN subpacking size M	15318392
PZN subpacking size L	15318386
PZN subpacking size XL	15318417
PZN subpacking size XXL	15318363
PZN subpacking size XXXL	-
measures & size	
length	400 mm
width	130 mm
heigth	45 mm
weights	

size	net weight	gross weight
XS	-	-
S	550 g	610 g
M	600 g	660 g
L	650 g	710 g
XL	700 g	760 g
XXL	750 g	810 g
XXXL	-	-

LOGISTIC DATA OUTER PACKING	
generell information	
material	carton
subpackings per outer packing	20
GTIN outer packing size XS	-
GTIN outer packing size S	4044941007298
GTIN outer packing size M	4044941007304
GTIN outer packing size L	4044941007311
GTIN outer packing size XL	4044941007410
GTIN outer packing size XXL	4044941014197
GTIN outer packing size XXXL	-
PZN outer packing size XS	-
PZN outer packing size S	-
PZN outer packing size M	-
PZN outer packing size L	-
PZN outer packing size XL	-
PZN outer packing size XXL	-
PZN outer packing size XXXL	-
measures & size	
length	460 mm
width	275 mm
heigth	410 mm
weights	

size	net weight	gross weight
XS	-	-
S	12.200 g	12.700 g
M	13.200 g	13.700 g
L	14.200 g	14.700 g
XL	15.200 g	15.700 g
XXL	16.200 g	16.700 g
XXXL	-	-

LOGISTIC DATA PALETTE	
general information	
kind of palett	euro-palette
measures & size	
cartons per layer	6
layers per palette	3
heigth of the palette	138 cm
weights	

size	net weight	gross weight
XS	-	-
S	229 g	254 g
M	247 g	272 g
L	265 g	290 g
XL	283 g	308 g
XXL	301 g	326 g
XXXL	-	-

AMPri Handelsgesellschaft mbH
Benzstraße 16
21423 Winsen (Luhe) - Germany
Telefon: +49 (0)4171 / 8480-0
Homepage: www.ampri.de
e-mail: info@ampri.de

Technical Data Sheet



Article-No.: **01240**

Description: **MED-COMFORT Blue Ultra 400**
Nitrile examination glove
blue, non sterile, powder free

WARNINGS AND SAFETY INFORMATION


storage / expiry date	Store gloves in original packaging in a cool and dry place without additional weight, protect from direct sunlight. Do not store near ozone sources (laser printers, copiers). The actual expiry time in use cannot be specified in general terms, as it depends on the general conditions of use. An individual risk assessment must be carried out in each case. The expiry date - valid for proper storage - is stated on the packaging.

use and control	Always use protective gloves only for the intended use and in the correct size. A check/risk assessment must be carried out to ensure that the gloves are suitable for the intended use, as the conditions at the workplace may deviate from those of the type test depending on temperature, abrasion and degradation. Breakthrough times and permeation levels are based on laboratory measurements and are determined using samples taken from the palm of the hand. The actual duration of protection of a glove with a specific substance can vary significantly due to the conditions of use (temperature, abrasion, stretching). In the case of aggressive chemicals, degradation (change in mechanical properties) can be an important factor to consider when selecting chemical-resistant gloves. This information does not reflect the actual duration of protection in the workplace and the distinction between mixtures and pure chemicals. The chemical resistance was determined under laboratory conditions only on the basis of samples from the palm and refers only to the chemicals tested. The situation may be different if the chemical is used in a mixture. The penetration resistance was evaluated under laboratory conditions and refers only to the tested specimen. The degradation results according to EN ISO 374-4 show the change in puncture resistance of the gloves after exposure to the tested chemical. Before use, the gloves must be checked for holes or damage.

disposal	Used gloves must be disposed of in accordance with the disposal regulations of the local waste disposal company. Unused gloves can be disposed of with household waste.

disinfection	Disinfection is not intended for these gloves and is the responsibility of the user.

warnings/ allergy information	Protective gloves are intended for single use only. This product contains dithiocarbamates, which may cause allergic reactions

donning and doffing instructions	

rev-no.: 8
date 03.11.2026

changes and errors excepted